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REPROCESSING INSTRUCTIONS

I. INTRODUCTION

The purpose of these reprocessing instructions is to stress the importance of proper handling, decontamination, cleaning, disinfection, sterilization and storage practices of surgical instruments. We would like to point out that the service life and functionality of a surgical instrument depend on the care given to all these steps.

An instrument that does not work properly is a potential risk for user and patient.

These instructions are for your information only. You have to comply with the statutory regulations in your country as well as the specific procedures that are in use in your hospital/clinic.

Our surgical quality stainless steel or titanium instruments meet very strict standards that have been defined by the European and international statutory authorities.

However, when exposed to inappropriate storage conditions, these instruments may corrode.

Therefore, a certain number of rules should be followed during their using cycle.

II. STORAGE CONDITIONS

The instruments should be stored:

- in appropriate, dry and clean rooms, without any temperature deviation,
- away from substances that may be corrosive and away from instruments with a corroded or damaged surface, preferably according to the items' type of material and without piling them up. Attention should be paid to their weight (put the heaviest instruments underneath the others).

III. CARE AND HANDLING

- Surgical instruments have a specific intended use (cutting, holding ...). Any improper use of the instrument will cause its deterioration.
- Before each use, make sure the instrument is clean and check its state under a microscope.
- All instruments must be wiped after each use, so as to eliminate any mark of antiseptic solution, any traces of tissues, cellular deposits, blood traces or any other deposits. Any substance that might prevent the decontamination liquid from entering hollow instruments must be removed as well.
- To return soiled instruments, please use appropriate trays or containers. The instruments should be returned in a liquid.
- Prolonged soaking in a physiological saline solution is prohibited. We recommend using distilled water.

IV. CLEANING INSTRUCTIONS

Validated decontamination and/or inactivation as well as cleaning procedures must be set up in collaboration with the pharmacist, the hygienist and the biomedical engineer of the various hospitals/clinics, in terms of efficacy, maintenance of the surgical instruments properties and compliance with national and international requirements (for France, please refer to circular letter DGS/5C/DHOS/E2/2001/138 dated 14 March 2001).

The above-mentioned procedures must be implemented immediately after using the instrument to avoid drying of any organic substance on the surface of the instrument or in its hollow parts.

According to the World Health Organization (WHO) classification, ophthalmic surgical procedure is performed on an organ with high infection risks due to the NCTA (non conventional transmissible agents). The instruments must therefore be handled separately.

The above-mentioned procedures should include the following conditions:

Immersion and decontamination

In order to prevent stains or impurities from drying and microbiological contamination on the instrument, it must be immersed immediately after use in an aldehyde-free decontaminating bath. This decontaminating immersion protects the staff in charge of instrument handling as well the environment and facilitates further cleaning procedures.

The following should then be taken into account:

- Comply with the instructions for use of the decontaminating solutions
- Minimum soaking time of 15 minutes
- Use a chloride-free solution
- The instruments must be immersed completely in a decontaminating bath, hinged-joint instruments must be completely opened and those made of several parts must be disassembled.
- The sieving trays in operating rooms that are exclusively dedicated to decontamination must be changed after each surgery.
- Rinse with distilled water.

Cleaning

Cleaning can be performed either in a machine or by ultrasound or manually if impurities are still visible on the instruments after decontamination.

Regardless of the cleaning method that is used, the detergents' instructions for use have to be followed as some of them can be used for decontamination and cleaning at the same time.

Manual cleaning

- Manual cleaning requires the use of a soft plastic brush for cleaning the handles and a soft clean cloth for the tips. Metal brushes are forbidden.
- Thoroughly rinse with distilled water

Mechanical cleaning

- The manufacturer's instructions are to be followed
- When positioning the instruments ensure that they do not come into contact with one another and that they are not damaged. Consider the weight and the shape of the instruments. Heavy instruments should be placed beneath lighter ones whereas hollow instruments must be put into a disinfection basket. Shadow areas must be avoided and hinged instruments must be completely opened and those made up of several parts must be disassembled.
- Do not put too many instruments in the disinfection baskets and in the machine.

Ultrasound cleaning

Ultrasound cleaning is not recommended.

However, if this cleaning method is applied, the following basic measures must be taken into account:

- The instructions given by the manufacturer must be followed.
- Temperature must be of 40-45°C.
- A frequency ranging from 25 to 50 KHz should be used.
- The ultrasound cleaning cycle should not exceed 5 minutes (excluding degassing time).
- Use a suitable cleaning and/or moisturising agent.
- Fragile instruments should be treated separately and the metal parts of the most delicate instruments should not be rubbed against metal.

NCTA inactivation

The agents that are responsible for Creutzfeldt-Jakob disease, non conventional transmissible agents (NCTA) also known as "prions", are particularly resistant to numerous physical and chemical treatments. Since the current treatments are not or only partially efficient (130°C/226°F temperature in humid environment, ultrasounds, UV, ethanol, formaldehyde ...), a particular inactivation procedure should be included in the pre-treatment procedure. For France, please refer to the recommendations provided by the French Ministry of Health concerning "the measures to be followed in the surgical field (...) with respect to the transmission risks of the Creutzfeldt-Jakob disease".

- Preferably use physical inactivation (between 134°C and 138°C for 18 minutes, autoclave sterilization for porous instruments)
- The use of sodium hypochlorite (bleach) is forbidden.
- In case of chemical inactivation, rinse with distilled water.

NB: for titanium instruments, treatment with sodium hydroxide may cause anodised blue to fade.

However it will not affect the state of the instrument's surface.

Despite chemical inactivation sterilization is always necessary.

Rinsing / Drying

All instruments must be rinsed with distilled water and dried rapidly and carefully with a clean cloth or with filtered compressed air.

V. CONTROL/ CLEANING/ REPAIRS

In order to prevent corrosion marks or stains from appearing, the instruments should be checked regularly and corrective measures should be taken during the instrument's life span. Check if the instruments are clean, in good condition and working properly.

Following control, some of the instruments may need lubrication, re-sharpening or repairing.

Lubrication can be performed in-house.

Yet, mind the quality of the lubricating solution which has to be water-permeable to permit autoclaving (in this case do not use silicone oil).

CROMA GmbH deals with the repairs of the instruments which it sells.

The instruments must be returned to the supplier after being cleaned according to all the above-mentioned instructions and in the same type of packaging as the one in which they were delivered in order to guarantee appropriate return conditions. Estimates are available on request. If your instrument cannot be repaired, we will contact you to discuss further steps.

VI. STERILIZATION

After placing the instruments in disinfection baskets (place the instruments in the first rank and make sure not to put too many), three sterilization methods may be used:

- Steam sterilization (porous load autoclave): 134 °C, for 18 minutes (theoretical time for ideal conditions)
- Steam sterilization (autoclave): 121 °C for 20 minutes.
- Ethylene oxide sterilization.

In case of steam sterilization, the instruments should not show any trace of oil or of any other impurity or contain any other correcting chemical solution.

Sterilization of an instrument in a closed position is forbidden.

The sterilization time and conditions recommended may vary from country to country.

Please refer to the regulations in your country and the specific procedures required by your hospital/clinic.

VII. RECOMMENDATIONS

- CROMA instruments that are delivered non-sterile must be decontaminated, cleaned, disinfected and/or sterilized before use. Once the instruments have been checked (no impurities visible), the first reprocessing cycle is used as a validation basis of the reprocessing procedure of all instruments.
- Rinsing the instruments with distilled water reduces the risks of stains, tarnishing and corrosion marks.

VIII. EXCHANGE/ GUARANTEE

- Only new instruments that have never been used and/or sterilized or those showing a defect when being used for the first time will be exchanged.
- All of our instruments have a one-year guarantee starting from the invoice date, provided that they are used and cleaned in accordance with the above-mentioned instructions and have not been dropped or subject to any other shock impact.

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